Endoform™ Antimicrobial

Restorative Bioscaffold

Dententor do Registro:

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Registro Anvisa: 80117580851

 **AROA BIOSURGERY LTD.**

2 Kingsford Smith Place
Airport Oaks, Auckland 2022
New Zealand

www.aroabio.com

AROA™, Endoform™ Antimicrobial and
AROA ECM™ are trademarks
of Aroa Biosurgery Ltd.

PI.9324.02 | ©March 2023

Instructions For Use

Endoform™ Antimicrobial

Restorative Bioscaffold

EN Product Description

Endoform™ Antimicrobial is a 100% w/w intact extracellular matrix (ECM), manufactured from ovine (sheep) forestomach tissue and incorporating approximately 12 ug/cm² (~0.3% w/w) ionic silver, intended to inhibit colonization of the device. This advanced wound care device retains the innate biological structure of the native tissue ECM, along with the associated macromolecules including elastin, fibronectin, glycosaminoglycans, and laminin. When rehydrated with exudate or sterile saline, **Endoform Antimicrobial** transforms into a soft conforming sheet, which is naturally incorporated into the wound over time.

Not made with natural rubber latex.

Prescription only.

Intended use

Endoform Antimicrobial is a sterile, single use ovine forestomach-derived extracellular matrix intended to cover, protect, and provide a moist wound environment.

Clinical benefit

Endoform Antimicrobial is intended for the management of acute and chronic wounds.

The extracellular matrix provides structure for the patient's own cells to build granulation tissue and supports epithelialization.

The device is effective against a broad spectrum of microbes, including *Acinetobacter baumannii*, *Candida albicans*, *Candida parapsilosis*, *Candida glabrata*, *Escherichia coli*, Methicillin resistant *Staphylococcus aureus* (MRSA), coagulase-negative *Staphylococci*, group A (beta-hemolytic) *Streptococci*, *Pseudomonas aeruginosa*, *Aspergillus niger*, and Vancomycin resistant *Enterococci* (VRE). The device provides sustained antimicrobial effectiveness within the dressing for up to 7 days.

Indications for use

Endoform Antimicrobial is indicated for the management of wounds including:

- Partial and full thickness wounds
- Pressure ulcers
- Venous ulcers
- Diabetic ulcers
- Chronic vascular ulcers
- Tunneled/ undermined wounds
- Surgical wounds (donor sites, grafts, post Moh's surgery, post laser surgery, podiatric, wound dehiscence)
- Traumatic wounds (abrasions, lacerations, second-degree burns, and skin tears)
- Draining wounds.

Contraindications

- **Endoform Antimicrobial** is derived from an ovine (sheep) source and should not be used on patients with known sensitivity to ovine (sheep) derived material or ionic silver.
- **Endoform Antimicrobial** is not indicated for use on third degree burns.

Risk statements

PRECAUTIONS

- **Endoform Antimicrobial** may be used when visible signs of infection are present in the wound area only when proper medical treatment addresses the underlying cause. **Endoform Antimicrobial** is not intended to be a substitute for appropriate treatment of infection.
- **Endoform Antimicrobial** is terminally sterilized via ethylene oxide. Do not use if the pouch is damaged.
- Discard the device if mishandling has caused possible damage or contamination.

- Single use product. Do not re-sterilize. Discard all unused portions. Device reuse may lead to wound infection.
- Always handle **Endoform Antimicrobial** using aseptic technique.
- Do not use past expiry date.
- **Endoform Antimicrobial** may cause transient discoloration of the wound bed and surrounding skin. Frequent or prolonged use of this product may result in permanent discoloration of the skin.
- Reduction in colonization or microbial growth on the device has not been shown to correlate with a reduction in infections in patients. Clinical studies to evaluate reduction in infection have not been performed.
- Healthcare Professionals should be aware that there is limited data on prolonged and repeated use of silver containing dressings, particularly in pediatrics and neonates.

Storage

Endoform Antimicrobial should be stored at room temperature in a clean and dry area, away from direct light.

Directions for use

These recommendations are designed to serve only as a general guideline. They are not intended to supersede institutional protocols or professional clinical judgment concerning patient care.

1. Wound Preparation

- a. Prepare the wound bed by cleansing, irrigation and if necessary, debride to remove debris, necrotic tissue or infected tissue.

2. Device Application

- a) Select a sheet of **Endoform Antimicrobial** which is slightly larger than the wound and apply aseptically.
- b) **Endoform Antimicrobial** can be applied as a whole sheet or trimmed so that it contacts the wound margins. Multiple sheets can be used to cover the entire wound bed.
- c) For ease of handling, apply the device by placing the dry material in the wound and rehydrating with exudate or sterile saline. When hydrated **Endoform Antimicrobial** transforms into a soft conforming sheet. Ensure that the device conforms to the underlying wound bed.
- d) To protect **Endoform Antimicrobial** from adhering to the secondary dressing, consider applying a non-adherent dressing over the device to help protect the tissue while facilitating an optimal moist wound healing environment.
- e) Secure the device using an appropriate secondary dressing. The secondary dressing should be changed according to standard of care and taking into account the level of exudate.
- f) **Endoform Antimicrobial** can be used in conjunction with compression therapy, off-loading, and negative pressure wound therapy.

3. Reapplication

- a) Duration of treatment and reapplication is determined by the physician and depends upon the wound type and conditions.
- b) The biodegradable **Endoform Antimicrobial** is absorbed and incorporated into the wound over time. It is recommended to assess the wound within 72 hours. Reapply when the device has been incorporated into the wound. **Endoform Antimicrobial** may last from 1 day and up to 7 days in the wound, depending on the wound type and environment.
- c) Carefully cleanse the wound surface in accordance with established procedures. Do not attempt to remove the off-white to golden gel that forms from residual **Endoform Antimicrobial** in the wound as it contains extracellular matrix components that assist in wound healing.
- d) If areas of the device are dry upon inspection, rehydrate with sterile saline and leave in place.
- e) It is not necessary to remove any residual **Endoform Antimicrobial** during dressing changes. However, if the product has been overlapped onto the peri-wound area, the remaining loose product that has not incorporated into the wound may be gently removed around the edges if desired.

Endoform™ Antimicrobial

Restorative Bioscaffold

- f) Change the secondary dressings as needed according to standard of care, taking into account the level of exudate, and when **Endoform Antimicrobial** is re-applied.

Appearance

Endoform Antimicrobial, like other silver-containing products, may darken upon storage, after hydration in saline, when exposed to light, or when in contact with body fluids and tissues. This darkening does not affect product performance.

Symbols Glossary

Symbols contained in ISO 15223-1 Medical devices – Symbols to be used with medical device labels, labeling and information to be supplied – Part 1: General requirements are indicated below.		
	5.1.1	Manufacturer
	5.1.3	Date of manufacture
	5.1.4	Use by/Expiration date
	5.1.5	Lot number
	5.1.6	Catalog number
	5.1.9	Importer
	5.2.3	Sterilized using ethylene oxide
	5.2.6	Do not resterilize
	5.2.8	Do not use if packaging is damaged
	5.2.11	Single sterile barrier system
	5.3.2	Keep away from sunlight
	5.3.4	Keep dry

	5.4.2	Do not reuse
	5.4.3	Consult instructions for use
	5.4.7	Contains a medicinal substance
	5.4.8	Contains biological material of animal origin
	5.7.7	Medical Device
	-	Number of units
	-	Fenestrated
	-	Non-Fenestrated
	-	High Flow
	-	Disc
Rx Only	-	Caution: Federal (USA) law restricts this device to sale by or on the order of a physician.

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