

SafeView

360° Panoramic Visualization Sterile Packaged • Fully Disposable





SafeView® Endoscopic Soft Tissue Release System



Next-Generation Features & Outcome-Driven Benefits

SafeView® Technology

- Transparent cannula
- Minimized cannula size
- Minimized incision

- 360° panoramic visualization
- Limited displacement to adjacent structures
- Easier, less disruptive insertion

Precision Control

- Independently operate arthroscope and knife blade within cannula
- Proprietary track technology
- Unlimited view on demand
- Precise and repeatable tissue release

Intuitive System

- Ergonomic instrument design
- Universal scope compatibility
- Simplified surgical steps
- Easily assimilated into any practice

Sterile Format

- Provided sterile, single-use
- Save time, save money, reduce infection potential

For more info please visit us at **SafeViewSurgery.com** or call **856.242.6979**

Carpal Tunnel Release

(1)

- Create a 1 cm transverse incision ulnar to the palmaris longus.
- Identify the volar forearm fascia, and incise in line with the skin incision.

 • Identify and protect the median nerve at this level.




- Insert the synovial elevator deep to the forearm fascia while remaining radial to the hook of the hamate.
 Palpate the distal edge of the transverse carpal
- A washboard effect will be felt as the elevator is moved longitudinally along the undersurface of




- Dilate the carpal tunnel space with the sequential dilators, aiming toward the third web space.

 Insertion depth is typically between 4–5 cm.









- Insert the SafeView® cannula, and palpate the palm to position the cannula just distal to the transverse
- Maintain posterior pressure on the hub of the cannula to preserve its position beneath the ligament.




- Insert a 4mm 30° standard arthroscope, and visualize the ligament and the deep fat distally as it overlaps the fibers of the transverse carpal ligament.


 Insert the In-situ rasp volar to the arthroscope with
- the handle oriented parallel to the forearm.



- Rotate 90 degrees to upright position.
- Rasp the undersurface of the transverse carpal ligament to clear away synovial tissue and improve visualization.








- Insert the forward cutting knife through the hub of the cannula.
- **(8)**
- Retract the skin proximally, and engage the proximal edge of the transverse carpal ligament.
- Divide the ligament under direct visualization.



- Following division, divide the volar antebrachial fascia with tenotomy scissors over a distance
- Skin closure is achieved in the usual fashion.
- Apply a soft dressing.







Cubital Tunnel Release

(1)

- Identify the medial epicondyle and olecranon process.
 A 2 cm longitudinal incision is made at the mid-point of these
- Divide the arcuate ligament (of Osborne), and identify the ulnar nerve
- A proximal release of the ulnar nerve may be completed under direct visualization.



 $(\mathbf{2})$

- To complete the distal release, first insert the synovial elevator deep to the forearm fascia.
- Insertion depth is typically between
 5-6 cm






 Dilate the cubital tunnel space with the sequential dilators.



(4)


- Insert the SafeView® cannula. Maintain posterior pressure on the hub of the cannula to preserve its position
- Insert a 4 mm 30 degree arthroscope. The cannula is appropriately positioned when the ulnar nerve can be visualized along the entire length of the cannula floor.






Optional:

- An optional rasp may be used to improve visualization of the forearm fascia prior to the release.
 Insert the rasp
- Insert the rasp superficial to the arthroscope with the handle oriented parallel to the elbow






 Rotate the rasp handle by 90 degrees to an upright position. Rasp the undersurface of the fascia to clear away any synovial tissue.



(7)

- Insert the forward cutting knife through the hub of the cannular
- Retract the skin proximally, and engage the proximal edge of the forearm fascia.
- Divide the fascia under direct visualization. Care must be taken to ensure that the ulnar nerve is visualized along the entire length of the cannula during division.




(8)

- Following a complete release of the cubital tunnel, the elbow is taken through a full range of motion. If ulnar nerve subluxation is detected, a transposition or epicondylectomy may be needed.
- Skin closure is achieved in the usual fashion.
- Per surgeon preference, apply a bulky soft bandage or long arm splint.



SafeView® Components



Ordering Information

Part Number	Description
1601.010	SafeView® Endoscopic Soft Tissue Release Kit

21st Century Solutions

Mission Surgical Innovations, LLC

992 Old Eagle School Road • Suite 907 • Wayne, PA 19087 O: 856.242.6979 • F: 856.295.8484 • E: sales@SafeViewSurgery.com

SafeViewSurgery.com